

## Photoelectric Effect

Rikuo Hasegawa  
Tutorial Group: C  
Lab Group: C2

February 6, 2019

Date of Experiment: February 6th, 2019

## 1 Introduction

In this experiment our goal is to measure the planck constant by examining the photoelectric effect using a halogen lamp and a PASCO h/e apparatus. We also want to measure the work function of the metal in the PASCO apparatus.

## 2 Theory

The energy of a photon is given in (eq. 1)

$$E = hf \quad (\text{eq. 1})$$

When photons with enough energy collide with a metal surface, it can cause electrons to be emitted. These electrons are sometimes referred to as photoelectrons. Each metal has a property which describes the amount of energy required before an electron gets emitted. This value is called the Work Function of the metal and is expressed using the symbol  $\Phi$ . The kinetic energy of a photoelectron is given in (eq. 2)

$$KE_{max} = hf - \Phi \quad (\text{eq. 2})$$

By measuring the voltage  $V_s$  of the emitted electrons, we can also find the kinetic energy as given in (eq. 3).

$$KE_{max} = eV_s \quad (\text{eq. 3})$$

where  $e$  is the charge of an electron, which is  $1.60217662 \times 10^{-19}$  [1].

## 3 Experimental Equipment and Method


### 3.1 Apparatus

Shown below is our experiment apparatus in Figure 1.

The equipment we used consists of several main components:

1. Halogen lamp as a source of white light
2. Convex lens to focus light onto PASCO apparatus
3. Filters of specific wavelength to limit frequency of light entering apparatus
4. PASCO h/e measurement apparatus
5. Digital Voltmeter to read values from PASCO apparatus

The uncertainty on the voltmeter is  $1 \times 10^{-3}$  [V]. We used filters of wavelength 400, 410, 440, 520, 546, and 600 [nm]. The uncertainty on these are  $5 \times 10^{-10}$  [m].



**Figure 1: Schematic of the apparatus used in the experiment**

Figure 1: [2]

### 3.2 Protocol

The experiment protocol is described below:

1. Turn on the halogen lamp and focus the light onto the PASCO device using the convex lens.
2. Turn on the PASCO device and the voltmeter.
3. Main loop, begin with 400 nm filter and move up:
  - (a) Insert the filter into the filter holder.
  - (b) Press the PRESS TO ZERO button on the PASCO device.
  - (c) Take measurement on voltmeter.
4. Repeat main loop using other filter wavelengths.
5. Turn all devices off.

## 4 Results


Our raw data is shown in Table 1. By plotting the  $eV_s$  against the frequency

wavelength [nm]	frequency [Hz]	Voltage [V]	$e \cdot V$ [J]
400	7.50E+14	1.87	3.00E-19
410	7.32E+14	1.91	3.06E-19
440	6.82E+14	1.69	2.70E-19
520	5.77E+14	1.50	2.40E-19
546	5.49E+14	1.32	2.12E-19
600	5.00E+14	0.90	1.44E-19

Table 1: Raw measurements

of light, we will be graphing (eq. 4), which can be derived from (eq. 2) and (eq. 3).

$$eV_s = hf - \Phi \quad (\text{eq. 4})$$

Graph 1:  $eV_s$  against frequency

The gradient of this graph should provide us with the Planck Constant and the y-intercept should provide us with the work function  $\Phi$  of the metal in the PASCO device.

## 5 Uncertainty Analysis

Using the regression feature of the data analysis toolkit in Microsoft Excel, we calculated the uncertainty in our measurements.

	<i>Value</i>	<i>Lower 95%</i>	<i>Upper 95%</i>
Work Function [J]	-1.10E-19	-2.68E-19	4.88E-20
Planck Constant [ $\text{m}^2\text{kg/s}$ ]	5.63E-34	3.14E-34	8.11E-34

Table 2: Uncertainty analysis

Thus, our measurements for the Planck Constant  $h$  and the work function  $\Phi$  are as follows:

$$h = 5.63 \pm 2.49 \times 10^{-34} [\text{m}^2 \text{kg/s}]$$

$$\Phi = -1.10 \pm 1.58 \times 10^{-19} [\text{J}]$$

The known value of the Planck Constant is  $h = 6.6260693 \times 10^{-34}$  [1]. Our measurements do contain this value, so we have successfully measured the

planck constant. However, with regard to both measurements, especially the work function, our random errors are very high.

## 6 Discussion and Conclusion

One of the most probable causes of the large uncertainty is that there were gaps in the filter for some of the filters, which would enable light of other frequencies to enter the PASCO device. This would render the measurements useless for the filters with this particular defect.

This experiment was successful in measuring the Planck Constant, albeit with a worryingly large uncertainty. With the work function of the metal, it is difficult to say that our experiment measured the value, as the work function cannot be negative.

## 7 Appendix

It is initially unintuitive that only frequency affects the voltage of photoelectrons. In order to rule out the effect of intensity of light, we took measurements with a filter which reduces the intensity of light regardless of wavelength. The observed voltage was not affected, as would be expected by the theory.

## References

- [1] P J Mohr. “The 2014 CODATA Recommended Values of the Fundamental Physical Constants”. In: (2016). URL: <http://physics.nist.gov/constants>.
- [2] UPCSE. *Physics Laboratory Experiments*. 2018, pp. 1–80.